


# Magic<sup>®</sup>



## FLOW-DEPENDENT MICROCATHETER SERIES

The **MAGIC** of true  
flow-dependent navigation



Now available through Balt USA

# The **ONLY** microcatheter available in a 1.2Fr OD for greater access through small distal vasculature



MODEL	REFERENCE	TOTAL LENGTH (cm)	GUIDEWIRE MAX (in)	DISTAL TIP OD (Fr)	MIN INTERNAL LUMEN (in)	OLIVE TIP MAX (in)	TURQUOISE PURSIL OD (Fr) x L (cm)	WHITE PURSIL OD (Fr) x L (cm)	PINK PURSIL OD (Fr) x L (cm)	GREEN TORVER OD (Fr) x L (cm)	DEAD SPACE (mL)
<b>MAGIC 1,2F</b>	MAGIC 1,2F	165	.008	1.2	0.009	-	1.2 x 12	1.8 x 3	2.4 x 30	2.7 x 114	0.38
<b>MAGIC 1,2F M</b>	MAGIC 1,2F M	165		1.2		-	1.2 x 3	1.8 x 12	2.4 x 30		0.38
<b>MAGIC 1,2F M 180</b>	MAGIC 1,2F M 180	180		1.2		-	1.2 x 3	1.8 x 17	2.4 x 40		0.40
<b>MAGIC 1,2F OLIVE</b>	MA 1,2F OLIVE	165		-		1.2	0.03	1.2 X 12	1.8 X 3	2.4 X 40	2.7 x 114
<b>MAGIC 1,5F</b>	MAGIC 1,5F	160	.009	1.5	0.011	-	-	1.5 x 15	2.4 x 25	2.7 x 114	0.37
<b>MAGIC 1,5F MP</b>	MA 1,5F MP	165		1.5		-	-	1.5 x 20			0.38
<b>MAGIC 1,5F OLIVE MP</b>	MA1,5F OLIVE MP	165		-		1.5	0.04	-	1.5 X 20		2.4 x 25
<b>MAGIC STD</b>	MAGIC STD	155	.010	1.8	0.013	-	-	1.8 x 10	2.4 x 25	2.7 x 114	0.37
<b>MAGIC MP</b>	MAGIC MP	165		1.8		-	-	1.8 x 20			0.38

MAGIC catheters are designed for general intravascular use. They may be used for the controlled, selective regional infusion of therapeutic agents or embolic materials into vessels.<sup>1</sup>

Federal (USA) law restricts this device to sale, distribution by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

1. Magic Catheters IFU – Ind 19

