

Hybrid

GUIDEWIRES

UNIQUE OPTIONS FOR COMPLEX ANATOMIES

Durability and Visibility


- Distal nitinol design maintains tip shape and enhances product durability
- Radiopaque distal tip ensures good visibility

Flexibility

- Available in straight or Double Angle shaped tips
- Broad options of ODs as small as .007"

Optimized Navigation

- Distal nitinol tip designed for shapeability and atraumatic navigation
- Proximal stainless steel for pushability and torqueability



ORDERING INFORMATION

Complete guidewire line in a stainless steel and nitinol design available in straight or double angle shape.

REFERENCE	DISTAL OD	PROXIMAL OD	LENGTH	DISTAL SHAPE
HYBRID007D	.007"	.012"	220cm	Straight
HYBRID007J	.007"	.012"	220cm	Double Angle
HYBRID008D	.008"	.012"	220cm	Straight
HYBRID008J	.008"	.012"	220cm	Double Angle
HYBRID10D300	.010"	.010"	310cm	Straight
HYBRID12D300	.012"	.012"	310cm	Straight
HYBRID1214D	.012"	.014"	200cm	Straight
HYBRID1214DA	.012"	.014"	200cm	Double Angle

Balt USA

29 Parker, Irvine, CA 92618 • tel 949.788.1443 fax 949.788.1444 • balt-usa.com

© 2018 BALT USA MKTG - 114 Rev. A

