

## GET THERE >>>

- Progressively softer coiled distal tip designed for safe and smooth navigation
- Large .088" ID to allow for delivery of a wide range of therapeutic devices

## STAY THERE >>>

- Braided reinforced proximal segment enhances support, pushability and kink resistance


# BALLAST

LONG SHEATH

**6F 088**

- Large .088” ID to allow for delivery of a wide range of therapeutic devices

**GET THERE. STAY THERE.**



	A	B	C	D	E
PART NUMBER	USABLE LENGTH	PROXIMAL OD	DISTAL OD	PROXIMAL ID	DISTAL ID
<b>BALLAST80</b>	80 cm	0.106"	0.100"	0.088"	0.088"
<b>BALLAST90</b>	90 cm				
<b>BALLAST100</b>	100 cm				

**Balt USA**

29 Parker, Irvine, CA 92618 • tel 949.788.1443 fax 949.788.1444 • [balt-usa.com](http://balt-usa.com)

© 2019 BAL T USA MKTG - 126 Rev. A

